

Formoterol

Catalog No: #S2580


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Formoterol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Adrenergic Receptor agonist
Calculated MW	344.4
Formulation	C ₁₉ H ₂₄ N ₂ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory system
 SMILES:C[C@H](CC1=CC=C(C=C1)OC)NC[C@@H](C2=CC(=C(C=C2)O)NC=O)O
 Pathways:GPCR/G
 ProteinReceptor:β2-adrenergic receptor
 Boiling pt:Melting pt:Solubility:DMSO: 10 mM
 Appearance:White Powder
 Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Kobayashi Y, et al. Pulm Pharmacol Ther. 2012 Jun;25(3):201-7.

Note: This product is for in vitro research use only