

Bedaquiline(TMC207; R207910)

Catalog No: #S2585


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bedaquiline(TMC207; R207910)
Brief Description	Inhibitors
Purification	98.00%
Target Name	Antibiotic inhibitor
Calculated MW	555.5
Formulation	C32H31BrN2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CN(C)CC[C@@](C1=CC=CC=C1)([C@H](C3=CC=CC=C3)C4=C(N=C5C=CC(=CC5=C4)Br)OC)O Pathways: Microbiology & virology Receptor: Tubercle bacillus Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Lakshmanan M, et al. J Young Pharm. 2013 Dec;5(4):112-115.

Note: This product is for in vitro research use only