

## Cabozantinib (XL184, BMS-907351)

Catalog No: #S2586


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cabozantinib (XL184, BMS-907351)
Brief Description	Inhibitors
Purification	99.71%
Target Name	c-Kit inhibitor; c-Met/HGFR inhibitor; TAM Receptor inhibitor; VEGFR inhibitor
Calculated MW	501.5
Formulation	C28H24FN3O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Cancer  
 SMILES: COC1=CC2=NC=CC(OC3=CC=C(C=C3)N(C(=O)C3(CC3)C(N)=O)C3=CC=C(F)C=C3)=C2C=C1OC  
 Pathways: Angiogenesis; Tyrosine Kinase/Adaptors/Receptor: Kit; c-Met; AXL; VEGFR2/KDR; VEGFR3/FLT4  
 Boiling pt: Melting pt: Solubility: DMSO: 100 mg/mL (199.39 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. You WK, et al. Y Res, 2011, 71(14), 4758-4768.

Note: This product is for in vitro research use only