

Ecabet Sodium

Catalog No: #S2587


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ecabet Sodium
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	402.48
Formulation	C ₂₀ H ₂₇ NaO ₅ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC(C)C1=C(C=C2C(=C1)CC[C@@H]3[C@@]2(CCC[C@@]3(C)C(=O)O)C)S(=O)(=O)[O-].[Na+] Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Lawrance IC. J Gastroenterol Hepatol. 2010 Jul;25(7):1182-4.

Note: This product is for in vitro research use only