

Veliparib (ABT-888)

Catalog No: #S2591


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Veliparib (ABT-888)
Brief Description	Inhibitors
Purification	98.92%
Target Name	PARP inhibitor
Calculated MW	244.29
Formulation	C13H16N4O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: C[C@@]1(CCCN1)C2=NC3=C(C=CC=C3N2)C(=O)N
 Pathways: Chromatin/Epigenetic; DNA Damage/DNA Repair
 Receptor: PARP1; PARP2
 Boiling pt: Melting pt: Solubility: DMSO: 17 mg/mL (69.58 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Albert JM, et al, Clin Y Res, 2007, 13(10), 3033-3042.

Note: This product is for in vitro research use only