

Eflornithine hydrochloride Monohydrate

Catalog No: #S2593


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Eflornithine hydrochloride Monohydrate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	236.64
Formulation	C ₆ H ₁₅ ClF ₂ N ₂ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: O.C1.NCCCC(N)(C(F)F)C(O)=O
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: >210°C
 Solubility: Soluble in DMSO
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Mohammed A, et al. Y Prev Res (Phila). 2014 Dec;7(12):1198-209.

Note: This product is for in vitro research use only