

Levomilnacipran

Catalog No: #S2594


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Levomilnacipran
Brief Description	Inhibitors
Purification	99.83%
Target Name	Norepinephrine inhibitor
Calculated MW	246.34
Formulation	C ₁₅ H ₂₂ N ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous

systemSMILES:CCN(CC)C(=O)[C@@]1(C[C@@H]1CN)C2=CC=CC=C2Pathways:NeuroscienceReceptor:5-HR-norepinephrine reuptakeBoling

pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Leo RJ, Brooks VL. Curr Opin Investig Drugs. 2006 Jul;7(7):637-42.

Note: This product is for in vitro research use only