

## GTBP Polyclonal Antibody

Catalog No: #40990

Package Size: #40990-1 50ul #40990-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GTBP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Mk
Specificity	GTBP Polyclonal Antibody detects endogenous levels of GTBP protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human GTBP.
Target Name	GTBP
Other Names	MSH6; GTBP; DNA mismatch repair protein Msh6; hMSH6; G/T mismatch-binding protein; GTBP; GTMBP; MutS-alpha 160 kDa subunit; p160
Accession No.	Swiss-Prot: P52701NCBI Gene ID: 2956
Uniprot	P52701
GeneID	2956;
SDS-PAGE MW	170kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of COS7 cells using GTBP Polyclonal Antibody



Western Blot analysis of A549 cells using GTBP Polyclonal Antibody

---

Note: This product is for in vitro research use only