

Thymopentin

Catalog No: #S2598


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thymopentin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	679.56
Formulation	C30H49N9O9
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:CC(C)[C@@H](C(=O)N[C@@H](CC=CC=C(C=C1)O)C(=O)O)NC(=O)[C@H](CC(=O)O)NC(=O)[C@H](CCCCN)NC(=O)[C@H](CCCN=C(N)N)NPathways:Endocrinology/HormonesReceptor:ERBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Zhu MX, et al. Exp Cell Res. 2015 Feb 15;331(2):387-98.

Note: This product is for in vitro research use only