

## Vidofludimus(4SC-101; SC12267)

Catalog No: #S2601


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Vidofludimus(4SC-101; SC12267)
Brief Description	Inhibitors
Purification	99.46%
Target Name	Dehydrogenase inhibitor
Calculated MW	355.36
Formulation	C20H18FNO4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism

systemSMILES:COC1=CC=CC(=C1)C1=CC(F)=C(NC(=O)C2=C(CCC2)C(O)=O)C=C1Pathways:MetabolismReceptor:Human DHODHBoling  
 pt:Melting pt:Solubility:Soluble in DMSOAppearance:Green PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Kulkarni OP, et al. Am J Pathol. 2010 Jun;176(6):2840-7.

Note: This product is for in vitro research use only