

Gα olf Polyclonal Antibody

Catalog No: #40991

Package Size: #40991-1 50ul #40991-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Gα olf Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	Gα olf Polyclonal Antibody detects endogenous levels of Gα olf protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Gα olf.
Target Name	Gα olf
Other Names	GNAL; Guanine nucleotide-binding protein G(olf) subunit alpha; Adenylate cyclase-stimulating G alpha protein; olfactory type
Accession No.	Swiss-Prot: P38405NCBI Gene ID: 2774
Uniprot	P38405
GeneID	2774;
SDS-PAGE MW	40kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of HEPG2-UV cells using Gα olf Polyclonal Antibody

Western Blot analysis of HepG2-UV cells using Ga α of Polyclonal Antibody



Note: This product is for in vitro research use only