


Description

Product Name	MSDC-0160
Brief Description	Inhibitors
Purification	99.53%
Target Name	Others
Calculated MW	370.42
Formulation	C19H18N2O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CCC1=CN=C(C=C1)C(=O)COC1=CC=C(CC2SC(=O)NC2=O)C=C1
 Pathways: Others
 Receptor: mTORT
 Boiling pt: Melting pt:
 Solubility: DMSO: 0.01 mg/mL (0.02 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
 Stock solution can be stored below -20 °C for several month.

References

1. Rohatgi N, et al. PLoS One. 2013, 8(5), e62012.

Note: This product is for in vitro research use only