

## Almorexant


Catalog No: #S2613

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Almorexant
Brief Description	Inhibitors
Purification	97.75%
Target Name	OX Receptor antagonist
Calculated MW	512.56
Formulation	C29H31F3N2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous

system SMILES: CNC(=O)[C@@H](C1=CC=CC=C1)N2CCC3=CC(=C(C=C3)C(F)(F)F)OC Pathways: GPCR/G

Protein; Neuroscience Receptor: OX1; OX2 Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Malherbe P, et al. Mol Pharmacol. 2009 Sep;76(3):618-31.

Note: This product is for in vitro research use only