

## Gα t1 Polyclonal Antibody

Catalog No: #40992

Package Size: #40992-1 50ul #40992-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Gα t1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Gα t1 Polyclonal Antibody detects endogenous levels of Gα t1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Gα t1.
Target Name	Gα t1
Other Names	GNAT1; GNATR; Guanine nucleotide-binding protein G(t) subunit alpha-1; Transducin alpha-1 chain
Accession No.	Swiss-Prot: P11488NCBI Gene ID: 2779
Uniprot	P11488
GeneID	2779;
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.

## Images



Western Blot analysis of COLO205 cells using Gα t1 Polyclonal Antibody

---

Note: This product is for in vitro research use only