

HBO1 Polyclonal Antibody

Catalog No: #40993

Package Size: #40993-1 50ul #40993-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	HBO1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF
Species Reactivity	Hu Ms Rt
Specificity	HBO1 Polyclonal Antibody detects endogenous levels of HBO1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human HBO1.
Target Name	HBO1
Other Names	KAT7; HBO1; HBOa; MYST2; Histone acetyltransferase KAT7; Histone acetyltransferase binding to ORC1; Lysine acetyltransferase 7; MOZ; YBF2/SAS3, SAS2 and TIP60 protein 2; MYST-2
Accession No.	Swiss-Prot: O95251NCBI Gene ID: 11143
Uniprot	O95251
GeneID	11143;
SDS-PAGE MW	75kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of Jurkat cells using HBO1 Polyclonal Antibody



Western Blot analysis of various cells using HBO1 Polyclonal Antibody

Note: This product is for in vitro research use only