

## AG-337;Nolatrexed dihydrochloride

Catalog No: #S2627


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	AG-337;Nolatrexed dihydrochloride
Brief Description	Inhibitors
Purification	99.50%
Target Name	DNA/RNA Synthesis inhibitor
Calculated MW	357.26
Formulation	C <sub>14</sub> H <sub>14</sub> Cl <sub>2</sub> N <sub>4</sub> O <sub>2</sub> S
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:CancerSMILES:CC1=C(SC2=CC=NC=C2)C2=C(NC(N)=NC2=O)C=C1Pathways:DNA Damage/DNA RepairReceptor:Thymidylate synthaseBoiling pt:Melting pt:301-302CSolubility:DMSO: > 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Hughes AN, et al. Br J Y. 2000 May;82(9):1519-27.

Note: This product is for in vitro research use only