

## Decernotinib(VX-509)

Catalog No: #S2636


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Decernotinib(VX-509)
Brief Description	Inhibitors
Purification	99.50%
Target Name	JAK inhibitor; Tyrosine Kinases inhibitor
Calculated MW	392.38
Formulation	C18H19F3N6O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:CancerSMILES:CC[C@](C)(C(=O)NCC(F)(F)F)NC1=NC(=NC=C1)C2=CNC3=C2C=CC=N3Pathways:Angiogenesis; Proteases/Proteasome; Chromatin/Epigenetic; JAK/STAT signaling; Stem CellsReceptor:JAK1; JAK2; JAK3; Tyk2Boling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Mahajan S, et al. J Pharmacol Exp Ther. 2015 Mar 11. pii: jpet.114.221176.

Note: This product is for in vitro research use only