

HDAC1 Polyclonal Antibody

Catalog No: #40994

Package Size: #40994-1 50ul #40994-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	HDAC1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	HDAC1 Polyclonal Antibody detects endogenous levels of HDAC1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HDAC1 around the non-phosphorylation site of S421.
Target Name	HDAC1
Other Names	HDAC1; RPD3L1; Histone deacetylase 1; HD1
Accession No.	Swiss-Prot: Q13547NCBI Gene ID: 3065
Uniprot	Q13547
GeneID	3065;
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using HDAC1 Polyclonal Antibody

Note: This product is for in vitro research use only