

Exherin (ADH-1)

Catalog No: #S2637


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Exherin (ADH-1)
Brief Description	Inhibitors
Purification	98.10%
Target Name	Others
Calculated MW	570.69
Formulation	C22H34N8O6S2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other
 SMILES: C[C@H]1C(=O)N[C@H](C(=O)N[C@@H](C(=O)N[C@H](C(=O)N1)CC2=CN=CN2)NC(=O)C(=O)N)C(C)C
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Perotti A, et al. Ann Oncol. 2009 Apr;20(4):741-5.

Note: This product is for in vitro research use only