


Description

Product Name	LY2811376
Brief Description	Inhibitors
Purification	99.77%
Target Name	BACE inhibitor; Beta Amyloid inhibitor
Calculated MW	320.36
Formulation	C15H14F2N4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:C[C@]1(CCSC(=N1)N)C2=C(C=C(C(=C2)C3=CN=CN=C3)F)FPathways:NeuroscienceReceptor:BACE1; AβBoling pt:Melting pt:Solubility:Ethanol: 64 mg/mL (199.77 mM); DMSO: 16 mg/mL (49.94 mM)Appearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. May P, et al. J Neurosci, 2011, 31(46), 16507-16516.

Note: This product is for in vitro research use only