

HDAC2 Polyclonal Antibody

Catalog No: #40995

Package Size: #40995-1 50ul #40995-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	HDAC2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	HDAC2 Polyclonal Antibody detects endogenous levels of HDAC2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HDAC2 around the non-phosphorylation site of S394.
Target Name	HDAC2
Other Names	HDAC2; Histone deacetylase 2; HD2
Accession No.	Swiss-Prot: Q92769NCBI Gene ID: 3066
Uniprot	Q92769
GeneID	3066;
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using HDAC2 Polyclonal Antibody



Western Blot analysis of Jurkat NIH-3T3 cells using HDAC2 Polyclonal Antibody

Note: This product is for in vitro research use only