

## INCB024360 analogue

Catalog No: #S2647


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	INCB024360 analogue
Brief Description	Inhibitors
Purification	98.71%
Target Name	IDO inhibitor
Calculated MW	271.64
Formulation	C9H7ClFN5O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism system  
 SMILES: C1=CC(=C(C=C1N/C(=C\2/C(=NON2)N)/N=O)Cl)F  
 Pathways: Metabolism  
 Receptor: IDO1  
 Boiling pt: Melting pt: Solubility: Ethanol: 54 mg/mL warmed (198.79 mM); DMSO: 54 mg/mL warmed (198.79 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Yue EW, et al. J Med Chem. 2009, 52(23), 7364-7367.

Note: This product is for in vitro research use only