


## Description

Product Name	GSK2656157
Brief Description	Inhibitors
Purification	99.50%
Target Name	ERK inhibitor
Calculated MW	416.45
Formulation	C23H21FN6O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:CancerSMILES:CN1C=C(C2=C1N=CN=C2N)C1=C(F)C2=C(C=C1)N(CC2)C(=O)CC1=CC=CC(C)=N1Pathways:MAPKReceptor:PERKBoling  
pt:Melting pt:Solubility:DMSO: 32 mg/mL (76.83 mM)Appearance:Off-white Color PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Atkins C, et al. Y Res, 2013, 73(6), 1993-2002.

Note: This product is for in vitro research use only