

Alda-1

Catalog No: #S2662


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alda-1
Brief Description	Inhibitors
Purification	99.50%
Target Name	Dehydrogenase agonist
Calculated MW	324.16
Formulation	C15H11Cl2NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
SMILES: C1C=CC(Cl)=C1C(=O)NCC1=CC=C(OCO2)C=C1
Pathways: Metabolism
Receptor: ALDH2
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several months.

References

1. Perez-Miller S, et al. Nat Struct Mol Biol. 2010 Feb;17(2):159-164.

Note: This product is for in vitro research use only