

Taladegib

Catalog No: #S2666


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Taladegib
Brief Description	Inhibitors
Purification	99.50%
Target Name	Hedgehog/Smoothened antagonist
Calculated MW	512.5
Formulation	C26H24F4N6O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:CancerSMILES:CN(C1CCN(CC1)C1=NN=C(C2=CC=NN2C)C2=CC=CC=C12)C(=O)C1=C(C=C(F)C=C1)C(F)(F)FPathways:GPCR/G Protein; Stem CellsReceptor:SmoothenedBoling pt:Melting pt:Solubility:Ethanol: 2 mg/mL (3.9 mM); DMSO: 0.66 mg/mL (1.28 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Redmond EM et al. Expert Opin Investig Drugs, 2011, 20(12),1649-1664.

Note: This product is for in vitro research use only