

## HDAC4 Polyclonal Antibody

Catalog No: #40997

Package Size: #40997-1 50ul #40997-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	HDAC4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	HDAC4 Polyclonal Antibody detects endogenous levels of HDAC4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HDAC4 around the non-phosphorylation site of S632.
Target Name	HDAC4
Other Names	HDAC4; KIAA0288; Histone deacetylase 4; HD4
Accession No.	Swiss-Prot: P56524NCBI Gene ID: 9759
Uniprot	P56524
GeneID	9759;
SDS-PAGE MW	119kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.


Not yet tested in other applications.

## Images



Western Blot analysis of HELA cells using HDAC4 Polyclonal Antibody

Western Blot analysis of HELA cells using HDAC4 Polyclonal Antibody



Note: This product is for in vitro research use only