

## Regadenoson

Catalog No: #S2671


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Regadenoson
Brief Description	Inhibitors
Purification	99.78%
Target Name	Adenosine Receptor agonist
Calculated MW	390.35
Formulation	C15H18N8O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Cardiovascular

SystemSMILES:CNC(=O)C1=CN(N=C1)C2=NC3=C(C(=N2)N)N=CN3[C@H]4[C@@H]([C@@H]([C@H](O4)CO)O)O

Pathways:GPCR/G  
ProteinReceptor:A2ABoling pt:Melting pt:158-160CSolubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Zoghbi GJ, Iskandrian AE. J Nucl Cardiol. 2012 Feb;19(1):126-41.

Note: This product is for in vitro research use only