


Description

Product Name	AS 602801
Brief Description	Inhibitors
Purification	99.92%
Target Name	JNK inhibitor
Calculated MW	457.55
Formulation	C25H23N5O2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: N#CC(C1=NC2=CC=CC=C2S1)C1=NC(OCC2=CC=C(CN3CCOCC3)C=C2)=NC=C1 Pathways: MAPK Receptor: JNK1; JNK2; JNK3
 Boling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Ferrandi C, et al. Mult Scler. 2011 Jan;17(1):43-56.

Note: This product is for in vitro research use only