

HDAC5/9 Polyclonal Antibody

Catalog No: #40998

Package Size: #40998-1 50ul #40998-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HDAC5/9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	HDAC5/9 Polyclonal Antibody detects endogenous levels of HDAC5/9 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HDAC5/9 around the non-phosphorylation site of S259/220.
Target Name	HDAC5/9
Other Names	HDAC5; KIAA0600; Histone deacetylase 5; HD5; Antigen NY-CO-9; HDAC9; HDAC7; HDAC7B; HDRP; KIAA0744; MITR; Histone deacetylase 9; HD9; Histone deacetylase 7B; HD7; HD7b; Histone deacetylase-related protein; MEF2-interacting transcription rep
Accession No.	Swiss-Prot: Q9UQL6/Q9UKV0NCBI Gene ID: 10014/9734
Uniprot	Q9UQL6
GeneID	10014;
SDS-PAGE MW	121kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using HDAC5/9 Polyclonal Antibody



Western Blot analysis of NIH-3T3 RAT-intestines cells using HDAC5/9 Polyclonal Antibody

Note: This product is for in vitro research use only