


## Description

Product Name	JNJ-1661010
Brief Description	Inhibitors
Purification	99.80%
Target Name	FAAH inhibitor
Calculated MW	365.45
Formulation	C19H19N5OS
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism

system SMILES: O=C(NC1=CC=CC=C1)N1CCN(CC1)C1=NC(=NS1)C1=CC=CC=C1 Pathways: Metabolism Receptor: FAAH (human); FAAH (rat) Boiling pt: Melting pt: Solubility: DMSO: 36 mg/mL (98.5 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Karbarz MJ, et al. Anesth Analg, 2009, 108(1), 316-329.

Note: This product is for in vitro research use only