


Description

Product Name	LY310762
Brief Description	Inhibitors
Purification	99.50%
Target Name	5-HT Receptor antagonist
Calculated MW	430.94
Formulation	C ₂₄ H ₂₇ FN ₂ O ₂ HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: Cl.CC1(C)C(=O)N(CCN2CCC(CC2)C(=O)C2=CC=C(F)C=C2)C2=CC=CC=C12
 Pathways: GPCR/G Protein;
 Neuroscience Receptor: 5-HT1D
 Boling pt: Melting pt: Solubility: DMSO: 0.4 mg/mL (0.92 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Pullar IA, et al. Eur J Pharmacol, 2004, 493(1-3), 85-93.

Note: This product is for in vitro research use only