

## Tranilast


Catalog No: #S2690

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tranilast
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others inhibitor
Calculated MW	327.33
Formulation	C18H17NO5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: COC1=CC(=C(C=C1)/C=C/C(=O)NC2=CC=CC=C2C(=O)O)OC  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: 166.2-168.2°C  
Solubility: DMSO: 66 mg/mL (201.63 mM)  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Suzawa H, et al. Jpn J Pharmacol, 1992, 60(2), 91-96.

Note: This product is for in vitro research use only