

## HDAC7 Polyclonal Antibody

Catalog No: #41000

Package Size: #41000-1 50ul #41000-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HDAC7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	HDAC7 Polyclonal Antibody detects endogenous levels of HDAC7 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HDAC7 around the non-phosphorylation site of S155.
Target Name	HDAC7
Other Names	HDAC7; HDAC7A; Histone deacetylase 7; HD7; Histone deacetylase 7A; HD7a
Accession No.	Swiss-Prot: Q8WUI4NCBI Gene ID: 51564
Uniprot	Q8WUI4
GeneID	51564;
SDS-PAGE MW	103kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of PC-3 cells using HDAC7 Polyclonal Antibody



Western Blot analysis of PC3 KB MCF7 HELA cells using HDAC7 Polyclonal Antibody

---

Note: This product is for in vitro research use only