

Asperosaponin

Catalog No: #S2698

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Asperosaponin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	929.1
Formulation	C47H76O18
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research

Area:CancerSMILES:C[C@@]12CCC3(CCC(CC3C1=CCC4C2(CCC5C4(CCC([C@@]5(C)CO)O[C@H]6[C@@H]([C@H]([C@H](CO6)O)O)C)C)C)C(=O)OC7[C@@H]([C@H]([C@@H]([C@H](O7)COC8[C@@H]([C@H]([C@@H]([C@H](O8)CO)O)O)O)O)O

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Song L, et al. Nat Prod Commun. 2014 Jun;9(6):773-8.

Note: This product is for in vitro research use only