

Mutant EGFR inhibitor

Catalog No: #S2705


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mutant EGFR inhibitor
Brief Description	Inhibitors
Purification	97.70%
Target Name	EGFR inhibitor
Calculated MW	520.03
Formulation	C27H30ClN7O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer SMILES: COC1=C(NC2=NC=C(Cl)C(=N2)C2=CNC3=CC=CC=C23)C=C(NC(=O)C=C)C(=C1)N(C)CCN(C)C Pathways: Angiogenesis; JAK/STAT signaling; Tyrosine Kinase/Adaptors Receptor: EGFR (mutant) Boling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only