

HDAC8 Polyclonal Antibody

Catalog No: #41001

Package Size: #41001-1 50ul #41001-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	HDAC8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	HDAC8 Polyclonal Antibody detects endogenous levels of HDAC8 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HDAC8 around the non-phosphorylation site of S39.
Target Name	HDAC8
Other Names	HDAC8; HDACL1; CDA07; Histone deacetylase 8; HD8
Accession No.	Swiss-Prot: Q9BY41NCBI Gene ID: 55869
Uniprot	Q9BY41
GeneID	55869;
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 293-UV cells using HDAC8 Polyclonal Antibody



Western Blot analysis of NIH-3T3 HeLa Jurkat cells using HDAC8 Polyclonal Antibody

Note: This product is for in vitro research use only