

D-Stachyose tetrahydrate

Catalog No: #S2713


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	D-Stachyose tetrahydrate
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	738.64
Formulation	C ₂₄ H ₅₀ O ₂₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other
 SMILES: O.O.O.O.OCC1OC(CO)(OC2OC(COC3OC(COC4OC(CO)C(O)C(O)C4O)C(O)C(O)C3O)C(O)C(O)C2O)C(O)C1O
 Pathways: Other
 sReceptor: Others
 Boling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: Faint Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Li W, et al. Mol Nutr Food Res. 2015 Nov 19.

Note: This product is for in vitro research use only