

1,3-Diphenylurea

Catalog No: #S2716


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	1,3-Diphenylurea
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	212.25
Formulation	C13H12N2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: O=C(NC1=CC=CC=C1)NC1=CC=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Huang W, et al. Bioorg Med Chem. 2010 Aug 1;18(15):5610-5.

Note: This product is for in vitro research use only