

## Salidroside

Catalog No: #S2717


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Salidroside
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	300.3
Formulation	C14H20O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: C1=CC(=CC=C1CCO[C@@H]2[C@@H]([C@H]([C@@H]([C@H](O2)CO)O)O)O)OP

Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Ethanol: 4 mg/mL (13.32 mM); Water: 60 mg/mL (199.8 mM); DMSO: 60 mg/mL (199.8 mM) Appearance: White

Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Hu X, et al. Cell Biol Toxicol, 2010, 26(6), 499-507.

Note: This product is for in vitro research use only