

Succinic acid


Catalog No: #S2719

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Succinic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	118.09
Formulation	C4H6O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:OC(=O)CCC(O)=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Sanjeeva Gandhi M, Mok YS. Chemosphere. 2014 Dec;117:440-6.

Note: This product is for in vitro research use only