

Tetramethylpyrazine

Catalog No: #S2722


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tetramethylpyrazine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	136.2
Formulation	C8H12N2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CC1=NC(C)=C(C)N=C1CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Guan D, et al. J Neurochem. 2015 Aug;134(3):551-65.

Note: This product is for in vitro research use only