

Ketoisophorone


Catalog No: #S2723

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ketoisophorone
Brief Description	Inhibitors
Purification	99.26%
Target Name	Others
Calculated MW	152.19
Formulation	C ₉ H ₁₂ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CC1=CC(=O)CC(C)(C)C1=O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Hegazy ME, et al. Z Naturforsch C. 2008 May-Jun;63(5-6):403-8.

Note: This product is for in vitro research use only