

Caryophyllene oxide

Catalog No: #S2724


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Caryophyllene oxide
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	220.36
Formulation	C15H24O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:C[C@@]12CC[C@@H]3[C@H](CC3(C)C)C(=C)CC[C@H]1O2Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Sain S, et al. Antiinflamm Antiallergy Agents Med Chem. 2014 Mar;13(1):45-55.

Note: This product is for in vitro research use only