

Sinomenine

Catalog No: #S2726


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sinomenine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	329.4
Formulation	C ₁₉ H ₂₃ NO ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CN1CC[C@@]23CC(=O)C(=C[C@@H]2[C@@H]1CC4=C3C(=C(C=C4)OC)O)O
 Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Ethanol: 8 mg/mL (24.28 mM); Water: 66 mg/mL (200.37 mM); DMSO: 66 mg/mL (200.37 mM)
 Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Liu L, et al. Int J Immunopharmacol, 1996, 18(10), 529-543.

Note: This product is for in vitro research use only