

Salvianolic acid B

Catalog No: #S2727


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Salvianolic acid B
Brief Description	Inhibitors
Purification	99.19%
Target Name	Others
Calculated MW	718.59
Formulation	C36H30O16
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: OtherSMILES: C1=CC(=C(C=C1C[C@H](C(=O)O)OC(=O)/C=C/C2=C3[C@H]([C@@H](OC3=C(C=C2)O)C4=CC(=C(C=C4)O)O)C(=O)O[C@H](CC5=CC(=C(C=C5)O)O)C(=O)O)O)OPathways: OthersReceptor: OthersBoiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Xue L, et al. Asian Pac J Trop Med. 2014 Apr;7(4):280-4.

Note: This product is for in vitro research use only