

Limonin

Catalog No: #S2728


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Limonin
Brief Description	Inhibitors
Purification	99.00%
Target Name	HIV Protease inhibitor
Calculated MW	470.52
Formulation	C26H30O8
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:C[C@@]12CC[C@H]3[C@]([C@@]14[C@H](O4)C(=O)O[C@H]2C5=COC=C5)(C(=O)C[C@@H]6[C@@]37COC(=O)C[C@@H]7OC6(C)C)CPATHWAYS:Microbiology&virology; Proteases/ProteasomeReceptor:HIV ProteaseBoiling pt:Melting pt:Solubility:DMSO: 44 mg/mL (93.51 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Battinelli L, et al. Planta Med, 2003, 69(10), 910-913.

Note: This product is for in vitro research use only