

## Scoparone

Catalog No: #S2729


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Scoparone
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	206.2
Formulation	C11H10O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:COC1=C(OC)C=C2C=CC(=O)OC2=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Lee SH, Jang HD. Exp Cell Res. 2015 Feb 15;331(2):267-77.

Note: This product is for in vitro research use only