

Usnic Acid

Catalog No: #S2731


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Usnic Acid
Brief Description	Inhibitors
Purification	99.98%
Target Name	Autophagy activator
Calculated MW	344.32
Formulation	C18H16O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:CC(=O)C1C(=O)C=C2OC3=C(C(C)=O)C(O)=C(C)C(O)=C3C2(C)C1=OPathways:AutophagyReceptor:AutophagyBoling
 pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it
 in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Su ZQ, et al. Int Immunopharmacol. 2014 Oct;22(2):371-8.

Note: This product is for in vitro research use only