

(-)-Epicatechin gallate

Catalog No: #S2732


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	(-)-Epicatechin gallate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	442.38
Formulation	C22H18O10
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:C1[C@H]([C@H](OC2=CC(=CC(=C21)O)O)C3=CC(=C(C=C3)O)O)OC(=O)C4=CC(=C(C(=C4)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 30 mg/mLAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Rosado H, et al. Int J Mol Sci. 2015 Jul 23;16(8):16710-27.

Note: This product is for in vitro research use only