

Sarsasapogenin

Catalog No: #S2733


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sarsasapogenin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Apoptosis inducer
Calculated MW	416.65
Formulation	C27H44O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:C[C@H]1CC[C@@]2([C@H]([C@H]3[C@@H](O2)C[C@@H]4[C@@]3(CC[C@H]5[C@H]4CC[C@H]6[C@@]5(CC[C@@H](C6)O)C)C)OC1Pathways:ApoptosisReceptor:Apoptosis inducerBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lim SM, et al. Int Immunopharmacol. 2015 Apr;25(2):493-503.

Note: This product is for in vitro research use only